
DRUG DELIVERY DEVICES

DELIVERING MORE.



Kimball Solutions produces the electronics for connected, high-tech pill bottles. These smart wireless pill bottles are currently being used by patients in pharmaceutical and research engagements. The bottles collect and send all adherence data in real-time. The system automatically analyzes this information and populates the data on the a secure dashboard. If doses are missed, patients can receive customizable alerts and interventions, using automated phone calls, text messages, and more. We take pride in producing assemblies for products that not only make the Medical industry better, but improve pharmaceutical research and patient care.

DRUG DELIVERY DEVICES

DELIVERING MORE.



We produce auto injectors and auto injector training devices for several medical customers.

With autoinjectors, the patient triggers the needle and the flow of medication through the application of pressure on the injection site.

Kimball Solutions:

- provides automated manufacturing services
- manufactures the electronics for the actuation of the drug delivery needle, corresponding sound prompts and connectivity
- manufactures the injection molded plastic components
- manages the cold chain process of pharmaceuticals
- completes the full assembly
- packages the devices for final distribution

Kimball Solutions continues to tailor our capabilities to meet the requirements and challenges of the medical market.